



Arterial catheter Kit

Precision with every beat. Confidence in every measurement.

Description

The MSB Arterial Catheter Kit was developed to ensure a continuous and accurate invasive reading of blood pressure and/or collection of samples for blood gas analysis, with total safety and reliability.

Ideal for procedures that require continuous hemodynamic monitoring, the catheter is inserted into a peripheral artery by percutaneous puncture and connected to a pressure transducer, coupled to a saline set.

The blood pressure oscillations captured at the catheter end are converted by the transducer into electrical signals, which are amplified and displayed on multiparameter monitors as curves and numerical values of arterial depression.

Specifications

- Guidewire: AISI 304 stainless steel
- Introducer needle: AISI 304 stainless steel
- Catheter: polyurethane, a flexible and biocompatible material that provides safety and performance during the procedure



Measurement specifications

Model	Device	Measurement
69005	KIT Arterial catheter	20Gx4,5CM FG 0,021x45CM AG 21GX4CM



+55 71 3024 4015
www.msbbrazil.com

Araponga Street, 364, Block 01, Lot 19,
Pitangueiras Lauro de Freitas - Bahia - Brazil
- 42701 330

✉ msb@msbbrazil.com
🌐 msb-medical-system-do-brasil



WARNING: Product mentioned in this material is intended for use by a physician. Before use, it is important to read the package insert thoroughly for the instructions for use, warnings, and possible complications associated with the use of this device.

Always check the regulatory status of the device in your region. It must not be reproduced, distributed, or extracted. Illustrations are artistic representations only and should not be considered as engineering drawings or photographs.

ARTERIAL CATHETER KIT - MSB - 80454419005 - 12/2025 - MSB - Version 1.2