



Ureteral Catheter

Double J



Description

The Double J MSB Ureteral Probe was developed with a focus on quality and well-being of the patient, it is used for urine drainage from the kidney and ureter. Made with high-tech soft polymer, it provides reduced pain, minimal trauma, and enhanced comfort during use. Its highly biocompatible material significantly reduces the risk of fouling, making the treatment safer and more effective.


Available in single packs or kit format with guidewire included, it adapts to the needs of the professional and the procedure. Always with a high standard of safety and efficiency.



Ureteral Catheter Double J

Specifications

- Customized measurements available on request;
- Manufactured with double sterile packaging, minimizing the risk of contamination;
- State-of-the-art technology to ensure optimal flexibility and best fit to the patient's urinary tract.



References and measures

Code	Diameter	Length	Kit	Guidewire
USDJ4712	4.7FR	12 cm	No	No
USDJ4714	4.7FR	12 cm	No	No
USDJ4716	4.7FR	16 cm	No	No
USDJ4720	4.7FR	20 cm	No	No
USDJ4722	4.7FR	22 cm	No	No
USDJ4724	4.7FR	24 cm	No	No
USDJ-4726	4.7FR	26 cm	No	No
USDJ-4728	4.7FR	28 cm	No	No
USDJ-472232	4.7FR	22-32 cm	No	No
USDJ-6026	6.0FR	26 cm	No	No
52038	4.7FR	26 cm	Yes	Teflon Straight 0.035"X150CM
USDJ-6026H35S	6.0FR	26 cm	Yes	Hydrophilic Straight 0.035"X150CM Standard
USDJ-472232H35S	4.7FR	22 - 32 cm	Yes	Hydrophilic Straight 0.035"X150CM
USDJ-602232H35S	4.7FR	22 - 32 cm	Yes	Hydrophilic Straight 0.035"X150CM
USDJ-4726H35S	4.7FR	26 cm	Yes	Hydrophilic Straight 0.035"X150CM Standard
USDJ-6026H35S	6.0FR	28 cm	Yes	Hydrophilic Straight 0.035"X150CM Standard



www.msbbrazil.com

Araponga Street, 364, Block 01, Lot 19,
Pitangueiras Lauro de Freitas - Bahia - Brazil
- 42701 330

 msb@msbbrazil.com

 msb-medical-system-do-brasil



WARNING: Product mentioned in this material is intended for use by a physician. Before use, it is important to read the package insert thoroughly for the instructions for use, warnings, and possible complications associated with the use of this device.

Always check the regulatory status of the device in your region. It must not be reproduced, distributed, or extracted. Illustrations are artistic representations only and should not be considered as engineering drawings or photographs.

Ureteral Probe - Double J - ANVISA Registration: 80454410004 - 12/2025 - MSB - Version 1.2