




# MSB Flex I

## Fiber Laser Cannula

- Laser amplifier connector compatibility

### Description

A reliable and safe solution for fragmentation of renal and ureteral stones. Through manufacturing with state-of-the-art technology, high-quality materials and rigorous validations, the result is the guarantee of an agile and successful surgical procedure.






# MSB Flex I Fiber Laser Cannula

## Characteristics

- Safe navigation: designed to provide safe navigation in the working channel during surgery;
- Excellent protection: it has additional protection that mitigates any energy leaks that can damage both the fiber and the laser;
- Wide compatibility: it has an SMA-905 connector, ensuring compatibility with most laser equipment available on the market.

## References

Code	Diameter	Length	Tip
62038	272 µm	3m	Flat
62040	365 µm	3m	Flat




[www.msbbrasil.com](http://www.msbbrasil.com)

Araponga Street, 364, Block 01, Lot 19,  
Pitangueiras Lauro de Freitas - Bahia - Brazil  
- 42701 330

✉ [msb@msbbrasil.com](mailto:msb@msbbrasil.com)

🌐 [msb-medical-system-do-brasil](http://msb-medical-system-do-brasil)



**WARNING:** Product mentioned in this material is intended for use by a physician. Before use, it is important to read the package insert thoroughly for the instructions for use, warnings, and possible complications associated with the use of this device.

Always check the regulatory status of the device in your region. It must not be reproduced, distributed, or extracted. Illustrations are artistic representations only and should not be considered as engineering drawings or photographs.

MSB Flex I - Laser Fiber Cannula - ANVISA Registration: 80454410012 - 12/2025 - MSB - Version 1.2