



# Guidewire MSB Hydro

## Description

The MSB HYDRO Guidewire is constructed with a single core that facilitates access to large stenotic lesions. Its extra lubricated hydrophilic coating offers greater agility and smoothness during procedures.



Single Core  
of Nitinol



## References


Code	Diameter	Length	Tip	Rigidity	Coating
UFGH-032150RHS	0.032"	150 cm	Straight	Standard	Hydrophilic
UFGH-035150RHS	0.035"	150 cm	Straight	Standard	Hydrophilic
UFGH-035150RHR	0.035"	150 cm	Straight	Rigid	Hydrophilic



[www.msbbrasil.com](http://www.msbbrasil.com)

Araponga Street, 364, Block 01, Lot 19,  
Pitangueiras Lauro de Freitas - Bahia - Brazil  
- 42701 330

 [msb@msbbrasil.com](mailto:msb@msbbrasil.com)  
 [msb-medical-system-do-brasil](http://msb-medical-system-do-brasil)



**WARNING:** Product mentioned in this material is intended for use by a physician. Before use, it is important to read the package insert thoroughly for the instructions for use, warnings, and possible complications associated with the use of this device.

Always check the regulatory status of the device in your region. It must not be reproduced, distributed, or extracted. Illustrations are artistic representations only and should not be considered as engineering drawings or photographs.

Hydrophilic Guidewire MSB - ANVISA Registration: 80454410001 - 12/2025 - MSB - Version 1.2