



MSB Catheter Introducer Set

Description

MSB Catheter Introducer Set is an indispensable device in any procedure. Its hemostatic valve technology allows introduction into both small and large vessels.

Specifications

- Ergonomically designed hub;
- Easy color identification on the hub and 3-way stopcock;
- Hydrophobic, kink-resistant, and flexible coating for smooth insertion;
- Atraumatic tip provides a smooth transition between the dilator and sheath during insertion;
- Locking dilator connector with French size indicator;
- Dilator with curved memory effect for easy sheath positioning;
- Allows for guidewires up to 0.035".



References

Code	Specifications
57109	Conj Intr Cateter MSB 12FR x 10cm D15cm FG 0.035" x 50 cm
57109	Conj Intr Cateter MSB 10FR x 10cm D15cm FG 0.035" x 50 cm
57304	Conj Intr Cateter MSB Magnum 18FR x 30cm D35cm



+55 71 3024 4015
www.msbbrazil.com

Araponga Street, 364, Block 01, Lot 19,
Pitangueiras Lauro de Freitas - Bahia - Brazil
- 42701 330

✉ msb@msbbrasil.com
🌐 msb-medical-system-do-brasil



WARNING: Product mentioned in this material is intended for use by a physician. Before use, it is important to read the package insert thoroughly for the instructions for use, warnings, and possible complications associated with the use of this device.

Always check the regulatory status of the device in your region. It must not be reproduced, distributed, or extracted. Illustrations are artistic representations only and should not be considered as engineering drawings or photographs.

MSB catheter introducer set - Anvisa Registration: 80454410009 - 12/2025 - MSB - Version 1.2