




M-Snare

Loop catheter



Pre-formed super elastic straps
constructed of nitinol and platinum.

Radiopaque marking at the distal tip.

Description

The M-Snare MSB catheter is designed for removal of foreign bodies within the peripheral arterial vasculature, providing precision and safety during the procedure. The bow is made of gold-plated tungsten, while the radiopaque marker of the catheter is made of platinum and iridium, providing clear visualization.

Measurement specifications

Model	Measurement	Loop
67009	6FRx100CM	15MMx120CM



+55 71 3024 4015
www.msbbrazil.com

Araponga Street, 364, Block 01, Lot 19,
Pitangueiras Lauro de Freitas - Bahia - Brazil
- 42701 330

✉ msb@msbbrazil.com
🌐 msb-medical-system-do-brasil



WARNING: Product mentioned in this material is intended for use by a physician. Before use, it is important to read the package insert thoroughly for the instructions for use, warnings, and possible complications associated with the use of this device.

Always check the regulatory status of the device in your region. It must not be reproduced, distributed, or extracted. Illustrations are artistic representations only and should not be considered as engineering drawings or photographs.

M-Snare - Anvisa Registration: 80454410015 - 12/2025 - MSB - Version 1.2