



# Centimeterized Angiographic Pigtail Catheter

## Description

The MSB Centimeterized Angiographic Catheter is an advanced instrument designed to deliver contrast media accurately and effectively at selected sites of the peripheral vascular system.

## Specifications

- Anatomical precision: the 20 platinum radiopaque markers, spaced every 1 cm, allow accurate anatomical measurements, contributing to detailed diagnoses;
- Controlled administration: the 10 lateral ports ensure precise and controlled delivery of contrast medium to the desired sites of the peripheral vascular system;
- Enhanced visibility;
- Versatility of use.

## Codes and references

Code	Measurement
66040	5FRX100CM 10F 20 B.RAD ANG CAT STRAIGHT



+55 71 3024 4015  
[www.msbbrazil.com](http://www.msbbrazil.com)

Araponga Street, 364, Block 01, Lot 19,  
Pitangueiras Lauro de Freitas - Bahia - Brazil  
- 42701 330

✉ [msb@msbbrazil.com](mailto:msb@msbbrazil.com)  
🌐 [msb-medical-system-do-brasil](http://msb-medical-system-do-brasil)



**WARNING:** Product mentioned in this material is intended for use by a physician. Before use, it is important to read the package insert thoroughly for the instructions for use, warnings, and possible complications associated with the use of this device.

Always check the regulatory status of the device in your region. It must not be reproduced, distributed, or extracted. Illustrations are artistic representations only and should not be considered as engineering drawings or photographs.

Centimeterized Pigtail - MSB - Anvisa Registration: 80454410014 - 12/2025 - MSB - Version 1.2