



UP!

Peripheral angioplasty balloon catheter


Description

The UP! Catheter Balloon is an advanced device developed to provide effective dilation in Transluminal-Percutaneous Angioplasty (PTA) procedures, with excellent performance in arterial stenosis.

Available in 0.018" and 0.035", the UP! Ensures optimal crossing, navigation and dilation, both in straight and tortuous vessels. Ideal for pre- and post-dilation of peripheral arteries.

Specifications

- Design developed for improved crossing and navigation
- Optimized tapered tip (0.040" in the 0.035" version)
- Compatible with 0.018" or 0.035" guidewire
- Distal hydrophilic coating facilitates vessel entry and exit





UP!

Peripheral angioplasty balloon catheter



Y-Manifold

The lateral lumen of the Y-Manifold is used to inflate and deflate the distal balloon.

Radiopaque Markers

Selective hydrophilic coating

Provides smooth delivery and pickup and enhanced traceability



Code	Diameter (mm)	Balloon length (mm)	Catheter length (cm)
0.018" guide wire			
73346	1,5	150	150
73339	2	120	150
73347	2	150	150
73363	2	200	150
73348	2,5	150	150
73341	3	120	150
73349	3	150	150
73365	3	200	150
73350	3,5	150	150
73343	4	120	150
73351	4	150	150
73367	4	200	150
73344	5	120	150
73352	5	150	150
73368	5	200	150

Code	Diameter (mm)	Balloon length (mm)	Catheter length (cm)
0.035" guide wire			
73160	5	40	135
73169	5	60	135
73192	5	120	135
73192	5	150	135
73207	5	200	135
73161	6	40	135
73170	6	60	135
73193	6	120	135
73198	6	150	135
73208	6	200	135
73162	7	40	135
73171	7	60	135
73194	7	120	135
73199	7	150	135
73209	7	200	135
73163	8	40	135
73172	8	60	135
73164	9	40	135
73173	9	60	135
73174	10	60	135
73166	12	40	135
73175	12	60	135



+55 71 3024 4015
www.msbbrazil.com

Araponga Street, 364, Block 01, Lot 19,
Pitangueiras Lauro de Freitas - Bahia - Brazil
- 42701 330

msb@msbbrazil.com

msb-medical-system-do-brasil



WARNING: Product mentioned in this material is intended for use by a physician. Before use, it is important to read the package insert thoroughly for the instructions for use, warnings, and possible complications associated with the use of this device.

Always check the regulatory status of the device in your region. It must not be reproduced, distributed, or extracted. Illustrations are artistic representations only and should not be considered as engineering drawings or photographs.

UP! Peripheral angioplasty balloon catheter - Anvisa Registration: 80454419008 - 12/2025 - MSB - Version 1.2