



MSB Ureteral

Introducer sheath

Description

The MSB ureteral introducer sheath has a hydrophilic body throughout its extension, ensuring excellent insertion and navigation. The conical dilator and soft polymer ensure a smooth advance in the ureter, with low risk of injury.

- HUB with snap-on
Easy and accurate




- Smooth transition between
the dilator and the sheath




MSB Ureteral Introducer sheath

Specifications

- Hydrophilic coating allows easy access and excellent handling;
- The cone-shaped dilator significantly improves access to the ureter, reducing the penetration effort and the risk of injury;
- Sheath with the atraumatic distal tip.



Hydrophilic
surface



Atraumatic tip

References

| Code | Diameter of the sheath | Length of the sheath |
|-------|---------------------------|-------------------------|
| 55004 | 10FR | 45 cm |
| 55005 | 10FR | 35 cm |
| 55011 | 10.7FR | 35 cm |
| 55012 | 10.7FR | 45 cm |



www.msbbrasil.com

Araponga Street, 364, Block 01, Lot 19,
Pitangueiras Lauro de Freitas - Bahia - Brazil
- 42701 330



msb@msbbrasil.com



msb-medical-system-do-brasil



WARNING: Product mentioned in this material is intended for use by a physician. Before use, it is important to read the package insert thoroughly for the instructions for use, warnings, and possible complications associated with the use of this device.

Always check the regulatory status of the device in your region. It must not be reproduced, distributed, or extracted. Illustrations are artistic representations only and should not be considered as engineering drawings or photographs.

Ureteral introducer sheath - ANVISA Registration: 80454419003 - 12/2025 - MSB - Version 1.2